

CERTIFICATE OF ANALYSIS

Product	accu-beads®+	Part Number	710258
----------------	--------------	--------------------	--------

Concentration Type	Acceptable Range	Lot #		Date of Manufacture	Expiration Date
35 M/ml	± 5.0 M/ml	241311351		2024.11.13	2026.11.13
Test Type and Quantity	20 µm Fixed Coverslip Chamber Calibration	S.D.	CV%	Calculated Hemacytometer Calibration Range*	
1 L - Bio	35.93	0.64	1.78	45.87-47.54	
1 L - QC	36.12	0.63	1.75	46.14-47.78	
5 ml - QC	35.46	0.82	2.3	45.03-47.16	
QC Analysis by: <i>Carlos A. Ulloa</i>		Reviewed by: Nicole Thibodeau			
Date: <i>2024-11-22</i>		Date: 2024-11-22			
Signature: <i>Carlos A. Ulloa</i>		Signature: Nicole Thibodeau		Digitally signed by Nicole Thibodeau Date: 2024.11.22 16:38:32 -05'00'	

* Calculation based on hemacytometer flow dynamics ratio as described in accu-beads Research Report 45. Please contact Hamilton Thorne for more information.