

**CERTIFICATE OF ANALYSIS**

<b>Product</b>	accu-beads®+	<b>Part Number</b>	710258
----------------	--------------	--------------------	--------

Concentration Type	Acceptable Range	Lot #		Date of Manufacture	Expiration Date
35 M/ml	± 5.0 M/ml	241405351		2024.05.14	2026.05.15
Test Type and Quantity	20 µm Fixed Coverslip Chamber Calibration	S.D.	CV%	Calculated Hemacytometer Calibration Range*	
1 L - Bio	35.48	0.36	1.02	45.65-46.59	
1 L - QC	34.82	0.76	2.19	44.27-46.26	
5 ml - QC	36.44	0.47	1.29	46.76-47.98	
QC Analysis by: <i>[Signature]</i>		Reviewed by: <i>Nicole Thibodeau</i>			
Date: <i>2024-06-24</i>		Date: <i>2024-06-25</i>			
Signature: <i>[Signature]</i>		Signature: <i>[Signature]</i>			

\* Calculation based on hemacytometer flow dynamics ratio as described in accu-beads Research Report 45. Please contact Hamilton Thorne for more information.