

CERTIFICATE OF ANALYSIS

Product Name	TRISTEL SPORICIDAL WIPE		
TRISTEL SPORICIDAL WIPE (Solution)			
Reference	TRL/018		
Date of Manufacture	2024/06/11		
Date of Expiry	2027/06/11		
This is to certify that Tristel Sporidical Wipe (Solution) bearing batch no: 816689 was tested on 2024/06/11 and found to comply with all requirements as set out in the specifications and method of analysis.			
PROPERTIES	SPECIFICATION	TEST METHOD	RESULTS
APPEARANCE	Clear and free of suspended matter	QTM029	Conforms
SPECIFIC GRAVITY	1.005 - 1.015	QTM006	1.006
pH	2.5 - 3.5	QTM007	2.9
TRISTEL SPORICIDAL WIPE (Sachet)			
Sachet batch numbers for this solution are as follows: 816689/1, 816689/2, 816689/3			

QUALITY CONTROL

Completed By		
Name	Signature	Date
ALEXANDRA FORTUNA	A. FORTUNA	12 Sep 2024
Reviewed By		
Name	Signature	Date
STEPH NAYLOR	[Signature]	12 Sep 2024

CERTIFICATE OF ANALYSIS

Product Name	TRISTEL ACTIVATOR FOAM		
TRISTEL ACTIVATOR FOAM (Solution)			
Reference	TRL/116 – FOAMER 007		
Date of Manufacture	2024/07/18		
Date of Expiry	2026/07/18		
<p>This is to certify that Tristel Activator Foam (Solution) bearing batch no: 816899 was tested on 2024/07/18 and found to comply with all requirements as set out in the specifications and method of analysis.</p>			
PROPERTIES	SPECIFICATION	TEST METHOD	RESULTS
APPEARANCE	Clear and free of suspended matter	QTM029	Conforms
SPECIFIC GRAVITY	1.000 – 1.010	QTM006	1.006
pH	10.0 – 11.5	QTM007	10.7
SODIUM CHLORITE	0.5 ± 0.05%	QTM004	0.52

QUALITY CONTROL

Completed By		
Name	Signature	Date
<i>ETOH LARA</i>	<i>[Signature]</i>	20 AUG 2024
Reviewed By		
Name	Signature	Date
Yason Osman	<i>[Signature]</i>	20 Aug 2024

QCA – Tristel Activator Foam- 10	Tristel Solutions Ltd.	Page 1 of 1
DRF7155		